



High Quality Imaging Solutions - American Service and Reliability

3016 NW 82nd Avenue Doral, Florida 33122 USA | +1 302 451.92.34 | info@medisono.com

www.medisono.com

Specifications subject to change without notice. For further details, please contact your MediSono sales representative.



P11 PLUS VET

Portable Colour Doppler Ultrasound System
The highest standard in shared- service ultrasound

3016 NW 82nd Avenue Doral, Florida 33122 USA | +1 302 451.92.34 | info@medisono.com

www.medisono.com

BRO-P11 PLUS Vet-Rev A-June-2021

P11 PLUS Vet

Medisono P11Plus Vet will provide the best clinical solution for veterinarians. P11Plus Vet is a high-end portable color Doppler system, it is incorporated with a lot of latest ultrasound technologies, such as Tissue Doppler Imaging, Anatomic M mode, and Panoramic Imaging and so on, which are rarely seen on veterinary ultrasound system.



The advanced imaging technologies, friendly ergonomic designs and powerful veterinary software packages make P11Plus Vet an unique veterinary software packages make P11VPlus Vet an unique professional VET ultrasound.

Medisono will commit to provide best products for the harmonious relationship between human beings and our animal friends.

Excellent Imaging and Processing Technologies

SNRT-scan Technology

With the probes to detect in more than one direction to eliminate artifacts, clinician can make clearer animal clinical diagnostic pictures.

Tissue Doppler Imaging

TDI is a new technique to display velocities of normal and pathologic myocardial structures during the cardiac cycle, facilitating clinical researchers to analyze and compare motions of different parts of animal hearts.

Steer M mode

By detecting ventricular wall motion in animals, clinicians and researchers can make more precise cardiac diagnosis and research reports.

Panoramic Imaging

Ultra-wide sight promises a clear vision and maximum of information on tissue and lesions, facilitating clinicians to make quick and precise diagnosis for large lesions.

The Main Features

- 15-inch anti-glare LCD monitor
- CFM, PDI, PW, CW, Steer M mode
- THI, TDI, Panoramic Imaging, Volumetric 4D Imaging
- IMT (intima-media thickness) measurement
- DICOM worklist, USB port, VGA

Professional Design for Veterinary Use

- Powerful veterinary measurement package
- User-definable veterinary exam modes
- Specialized veterinary body marks
- Professional veterinary transducers

Other Strength Features

- Two standard transducer sockets compatible with most color Doppler transducers of Medisono
- Built-in high capacity transducers of Medisono
- High mobility, slim and compact trolley
- Abundant and professional veterinary transducers

